



NYUCancerInstitute

INSTRUCTIONS FOR USE OF CADD INFUSION PUMP

1) CHANGE YOUR BATTERY:

Change your battery only if you hear beeping, and the screen says "LO BAT". Pump uses two AA batteries. Remove cover over battery chamber, located on the back of the pump, by sliding it off. Replace batteries following the illustration inside the empty battery chamber. Replace cover. **Remember to restart pump as below.**

2) STARTING AND STOPPING PUMP:

You may need to stop the pump, only if directed to do so by medical personnel. To stop the pump press down on "START/STOP" key and hold down until three dashes(---) appear on screen. Be sure to apply constant pressure until all three dashes appear. The word "STOPPED" will appear on the screen. To restart your pump (i.e. after changing battery or pump becomes turned off inadvertently) press down on "START/STOP" key and hold down until the three dashes(---) disappear from screen. When pump is running, the words "RUN/RES VOL" will be on the screen.

3) ALARMS AND TROUBLE SHOOTING:

When the reservoir volume is low (less than 5ml) the pump will emit a short beep every four to five minutes. This means your cassette should be changed very soon. This should only occur the day you are scheduled for a cassette change. If this happens when you don't have an appointment to change the cassette, please call your doctor.

If the battery is low, you will hear three beeps every 5 minutes. Change the battery as directed above.

Sometimes tubing may become kinked or one of the three small white plastic tubing clamps may be in the closed position. In this case, you will hear a constant tone and the screen will display "HIGH PRESSURE". First check all the clamps to be sure they are not pinching the tubing. If so, unclamp the tubing. **Then remember to restart your pump.**